

Homogenizers



X 1000 Unidrive
X 1000D Unidrive



X 120



X 1740



Please note that the maximum
depth of immersion is:

fine **F**

low viscosity **N**

viscous **V**

knife **M**


T 6	T 6 long	T 10	T 17	T 20	G 20	T 30	G 30	T 40	G 40
2,0 - 50 ml		10,0 - 250 ml	250 - 1500 ml	250 - 5000 ml		500 - 5000 ml		5000 - 20000 ml	
55 mm	120 mm	150 mm	150 mm	200 mm	215 mm	250 mm	280 mm	250 mm	280 mm
60420-000F	60410-000F	60421-000F	60426-000F	60422-000F	60423-000F	60427-000F	60428-000F	60414-000F	60415-000F
60420-000N	60410-000N	60421-000N	60426-000N	60422-000N	60423-000N	60427-000N	60428-000N	60414-000N	60415-000N
60420-000V	60410-000V	60421-000V	60426-000V	60422-000V	60423-000V	60427-000V	60428-000V	60414-000V	60415-000V
60420-000M	60410-000M	60421-000M	60426-000M	60422-000M	60423-000M	60427-000M	60428-000M	60414-000M	60415-000M

T: стандартное применение
(например, для жидкостей)

G: для агрессивных сред или
если требуется давление




Overhead Stirrer Accessories





Part No.

60617-0050	Dissolver ø 100 mm + shaft
60617-0000	Dissolver ø 100 mm
60617-0010	Dissolver shaft ø 36 mm
60618-0000	Dissolver shaft ø 8x350 mm
60617-0060	Dissolver ø 36 mm + shaft



Part No.

60618-0050	Stirring blade + shaft
60619-0000	Stirring blade (L = 90 mm)
60619-0010	Stirring blade (L = 150 mm)
60618-0000	Stirring shaft (ø 8x350 mm)
60615-0040	Stirring blade tool horizontal (ø 4x230 mm)
60615-0020	Stirring blade tool vertical (ø 4x250 mm)



Part No.

Blade Stirrer	Length including blade	Diameter of rod	Diameter of blade
8B00940400	400 mm	8 mm	70 mm
8B00940800	500 mm	10 mm	70 mm
8B00950400	400 mm	8 mm	70 mm
8B00950800	500 mm	10 mm	70 mm

Part No.

Centrifugal Stirrer	Length including blade	Diameter of rod	Diameter of blade
8B00921000	400 mm	8 mm	90/15 mm

Part No.

Blade Stirrer	Length including blade	Diameter of rod	Diameter of blade
8B00970000	300 mm	6 mm	50 mm
8B00970300	300 mm	8 mm	50 mm
8B00970400	400 mm	8 mm	50 mm
8B00970500	500 mm	8 mm	50 mm

Part No.

Propeller Stirrer	Length including blade	Diameter of rod	Diameter of blade
8B00910000	300 mm	6 mm	50 mm
8B00910300	300 mm	8 mm	50 mm
8B00910400	400 mm	8 mm	50 mm
8B00910500	500 mm	8 mm	50 mm
8B00911000	400 mm	8 mm	70 mm
8B00911400	500 mm	10 mm	100 mm
8B00900000	300 mm	6 mm	50 mm
8B00900400	400 mm	8 mm	50 mm
8B00901000	400 mm	8 mm	100 mm